

Anti-Tipper

Wheel Type	Seat to Floor Height	Anti-Tipper Part Number (1 each)	Assembly Quantity Per Wheelchair
12COM	12-1/2" - 13-1/2"	1159508	2
12COM	15-1/2" - 16-1/2"	1163441	2
12COM	17" - 19"	1163442	2
16COM	12-1/2" - 15"	1159508	2
16COM	17" - 19"	1163442	2
18COM	13-1/2" - 15"	1159508	2
18COM	15-1/2" - 16"	1163441	2
18COM	18" - 19"	1163442	2
20COM	ALL	1163441	2
22COM and 24COM	ALL	1163442	2