







Bodypoint Side-Release Ankle Huggers and Hardware




Installation Instructions



(4) Hard Plastic Loops



(4) Adjustment Buckles



(8) Screws

Bodypoint Side-Release Ankle Huggers and Hardware

Description	Part Number	Circumference Above Ankle Adjustment Range
Kit, Side-Release Ankle Huggers and Hardware, Small	1163628 <small>See Note(s) A</small>	6-1/2" - 8"
Kit, Side-Release Ankle Huggers and Hardware, Medium	1163629 <small>See Note(s) A</small>	7-1/2" - 9"
Kit, Side-Release Ankle Huggers and Hardware, Large	1163630 <small>See Note(s) A</small>	8-1/2" - 11-1/2"
Kit, Side-Release Ankle Huggers and Hardware, Extra Large	1163631 <small>See Note(s) A</small>	11" - 13"
NOTE: A - Kit includes the right and left Side Release Ankle Huggers, attaching hardware and installation instruction sheet. Individual components are non-serviceable.		