




# Canes

Quad - Model 7823, 7827, 7829, 3920-2

Model 7823, Medium Base  
Model 7829, Large Base



Model 3920- 2, Medium Base  
Model 7827, Large Base



# Canes

## Quad - Model 7823, 7827, 7829, 3920-2

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B,H	1048091-NLA	NLA - Hand Grip, Standard	1	1
1	A,B	1130776-NLA	NLA - Kit, Standard Grips	10	1
1	C,D	1062839-NLA	NLA - Hand Grip, Dual Cam	1	1
1	C,D,I	1130778-NLA	NLA - Kit, Dual Cam Grips	10	1
2		1053041-NLA	NLA - Spring, Torsion w/ Double Buttons	1	1
3		1055948	Kit, Camlock	1	1
			Lever, Cam, Black	1	1
			Housing, Cam	1	1
			Instruction Sheet, Camlock	1	1
4	A,C,E	1052368-NLA	NLA - Base Assembly, Narrow	1	1
4	B,D,E	1052369-NLA	NLA - Base Assembly, Large	1	1
5	A,C,F,J	1049993-NLA	NLA - Tip, Cane, Rubber, Black (1/2" Tubing)	1	4
5	B,D,G	1049992-NLA	NLA - Tip, Cane, Rubber, Black (5/8" Tubing)	1	4
<p>NOTE: A - To be used on Model 3920.            B - To be used on Model 7827.            C - To be used on Model 7823.            D - To be used on Model 7829.            E - Individual item no longer available, order complete cane for replacement.            F - For 1/2" Tubing.            G - For 5/8" Tubing.            H - Individual item no longer available, refer to Kit Part #1130776.            I - Individual item no longer available, refer to Kit Part #1130778.            J - After 8/28/08 appearance of the rubber tip changed. New tips are functionally equivalent, but should be ordered in quantities of 4 if uniform appearance is required or desired.</p>					