



MKIV Controller, Programmer and Mounting Hardware

NOTE: See the MKIII Controller Mounting Hardware page for information on mounting the MKIV controller on the back canes.



Battery Box Top Mount



Back Cane Mount

MKIV Controller, Programmer and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1039860	Screw, Hex Head w/ Patch (5/16-18 x 3/4")	1	2
2		1037943-NLA	NLA - Bracket with Keepers, Controller Mounting	1	1
3		1026158	Locknut (5/16-18)	1	2
4		1021260	Tie Wrap (5-5/8" L)	1	1
5		1105707-NLA	NLA - Controller, MKIV Quad Pulse	1	1
5	A,G	1095269-NLA	NLA - Controller, MKIV RII, Retro	1	1
6		16348X004-NLA	NLA - Rivet	1	4
7	H	50000X025-NLA	NLA - Clip, Lock Over Latch	1	2
8		1013149-NLA	NLA - Label, Caution	1	1
9		1027239-NLA	NLA - Bracket, Controller Mounting, Black	1	1
10	B	1109091-NLA	NLA - Programmer, Remote	1	1
11		1020818	Switch, Red Flex Ribbon (Tape) Reset	1	1
11	I	1123276	Kit, Switch, Black Egg Reset w/ Hook/Loop Adhesive Back Fastener & Attaching Screws (EGSBLK)	1	1
	C,D	1064404-NLA	NLA - Virtual Service Modem Kit	1	1
	C,E	1091631-NLA	NLA - Virtual Service Software Kit	1	1
	C,F	1046117-NLA	NLA - Printer and Interface	1	1
	C	1095204-NLA	NLA - Cable, Virtual Service, Charger/Modem	1	1
	C	1091980-NLA	NLA - Disk, 3.5 H.D. Floppy with V.S. Software	1	1
	C	1091637-NLA	NLA - Adaptor, Null Modem	1	1
	C	1091636-NLA	NLA - Service Manual - V.S. Software Kit	1	1
	C	1026572-NLA	NLA - Paper, Printer	1	1
	C	1027280-NLA	NLA - Harness, Printer Wiring	1	1
	C	1092610-NLA	NLA - Joystick, RII	1	1

NOTE: A - RII retro controller after 4/30/00 and RII joystick part number 1092610 are direct replacements for joystick/controller part number 1080723.
 B - Remote Programmer part number 1109091 with the connectors (3-pin and 5-pin) is the only programmer functionable with the new RII joystick part number 1092610.
 C - Not Shown
 D - Part number 1095204 is available for individual sale.
 E - Part numbers 1095204, 1091980, 1091637, and 1091636 are available for individual sale.
 F - Part numbers 1026572 and 1027280 are available for individual sale.
 G - After 4/30/00
 H - For Keeper
 I - Not Shown. After 5/23/04 the Egg Reset Switch became available as a service part and can be used as a direct replacement for the tape reset switch part # 1020818.