


Canes

Composite Quad Cane

Model 7821 - Adult



Canes

Model 7821 - Adult


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1048091-NLA	NLA - Hand Grip, Standard	1	1
1		1130776-NLA	NLA - Kit, Standard Grips	10	1
2		1053041-NLA	NLA - Spring, Torsion w/ Double Buttons	1	1
3		1055948	Kit, Camlock	1	1
			Lever, Cam, Black	1	1
			Housing, Cam	1	1
			Instruction Sheet, Camlock	1	1
4	B	1049993-NLA	NLA - Tip, Cane, Rubber, Black (1/2" Tubing)	1	5
5		18860X007	Screw, Phillips Pan Head (1/4-20 x 2-1/4")	1	1

NOTE: A - Individual item no longer available, refer to Kit Part #1130776.
 B - After 8/28/08 appearance of the rubber tip changed. New tips are functionally equivalent, but should be ordered in quantities of 5 if uniform appearance is required or desired.

Canes


Orthopedic Grip and Folding Cane

Model 8928 or 8929



Single Side Shown

Model 40800-6



Canes

Model 8928, 8929 and Model 40800-6

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1052204-NLA	NLA - Clip, Spring	1	1
2	A,B,D	1049735-NLA	NLA - Tip, Rubber, Black (3/4" Tubing)	1	1
2	C,E	1065416	Tip, Rubber, Black (3/4" Tubing)	1	1
NOTE: A - To be used with Model 8928 and 8929. B - Due to style changes, when replacing rubber tip, measure inside diameter. C - To be used with Model 40800. D - Item no longer available, refer to part #1048113 for replacement. E - After 10/1/08, appearance of the rubber tip has changed. New tips are functionally equivalent.					