

9000XDT & 9000XT Recliner Conventional Adjustable Height Arms (19JY and 29JY)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	8881130809	Kit, Armrest, Full Length, Black - Right (Style 19JY)	1	1
2			Spring, Torsion with Button Lock	1	2
3			Plug, Armrest, Black (3/4") (18GA)	1	2
4			Rivet	1	6
5			Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
6			Pin, Locking	1	1
7			Spring, Arm Lock	1	1
8			Handle, Lock, Black - Right	1	1
9			Rivet, Tubular	1	1
10			Guard, Clothing	1	1
11			Arm Assembly, Adjustable Height (2" Lower)	1	1
12			Anti-Rattle	1	2
13			Handle, Quad Release	1	2
14			Tube, Armrest, Full Length, Black (2" Lower)	1	1
15			Plug Button, Black (7/8")	1	2
17			NLA - Pad, Arm, Full Length w/ Law Label	1	1
16			Plug Button, Black (5/8")	1	2
18			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	1
4		1028645	Rivet	1	6
10		1068511	Guard, Clothing	1	1
1	A,D	8881130810	Kit, Armrest, Full Length, Black - Left (Style 19JY)	1	1
2			Spring, Torsion with Button Lock	1	2
3			Plug, Armrest, Black (3/4") (18GA)	1	2
4			Rivet	1	6
5			Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
6			Pin, Locking	1	1
7			Spring, Arm Lock	1	1
8			Handle, Lock, Black - Left	1	1
9			Rivet, Tubular	1	1
10			Guard, Clothing	1	1
11			Arm Assembly, Adjustable Height (2" Lower)	1	1
12			Anti-Rattle	1	2
13			Handle, Quad Release	1	2
14			Tube, Armrest, Full Length, Black (2" Lower)	1	1
15			Plug Button, Black (7/8")	1	2
16			Plug Button, Black (5/8")	1	2
17			NLA - Pad, Arm, Full Length w/ Law Label	1	1
18			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	1
4		1028645	Rivet	1	6
10		1068511	Guard, Clothing	1	1
1	A,D	8881130807	Kit, Armrest, Desk Length, Black - Right (Style 29JY)	1	1
2			Spring, Torsion with Button Lock	1	2
3			Plug, Armrest, Black (3/4") (18GA)	1	2

9000XDT & 9000XT Recliner Conventional Adjustable Height Arms (19JY and 29JY)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
4			Rivet	1	6
5			Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
6			Pin, Locking	1	1
7			Spring, Arm Lock	1	1
8			Handle, Lock, Black - Right	1	1
9			Rivet, Tubular	1	1
10			Guard, Clothing	1	1
11			Arm Assembly, Adjustable Height (2" Lower)	1	1
12			Anti-Rattle	1	1
13			Handle, Quad Release	1	1
14			Tube, Armrest, Desk Length, Black (2" Lower)	1	1
15			Plug Button, Black (7/8")	1	2
16			Plug Button, Black (5/8")	1	2
17			Arm Pad, Desk Length w/ Law Label	1	1
18			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
4		1028645	Rivet	1	6
10		1068511	Guard, Clothing	1	1
1	A,D	8881130808	Kit, Armrest, Desk Length, Black - Left (Style 29JY)	1	1
2			Spring, Torsion with Button Lock	1	2
3			Plug, Armrest, Black (3/4") (18GA)	1	2
4			Rivet	1	6
5			Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
6			Pin, Locking	1	1
7			Spring, Arm Lock	1	1
8			Handle, Lock, Black - Left	1	1
9			Rivet, Tubular	1	1
10			Guard, Clothing	1	1
11			Arm Assembly, Adjustable Height (2" Lower)	1	1
12			Anti-Rattle	1	1
13			Handle, Quad Release	1	1
14			Tube, Armrest, Desk Length, Black (2" Lower)	1	1
15			Plug Button, Black (7/8")	1	2
16			Plug Button, Black (5/8")	1	2
17			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
18			Arm Pad, Desk Length w/ Law Label	1	1
4		1028645	Rivet	1	6
10		1068511	Guard, Clothing	1	1
19	B,D	8881127509	Kit, Arm Pad w/Screws - Full Length	1	1
17			NLA - Pad, Arm, Full Length w/ Law Label	1	1
18			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
19	B,D	8881127510	Kit, Arm Pad w/Screws - Desk Length	1	1
17			Arm Pad, Desk Length w/ Law Label	1	1
18			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
	C	7100053	Arm Lock Assembly, Black - Right	1	1

9000XDT & 9000XT Recliner Conventional Adjustable Height Arms (19JY and 29JY)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	C	7100054	Handle, Lock, Black - Right	1	1
			Rivet, Tubular	1	2
			Pin, Locking	1	2
			Spring, Arm Lock	1	2
			Arm Lock Assembly, Black - Left	1	1
			Handle, Lock, Black - Left	1	1
			Rivet, Tubular	1	2
			Pin, Locking	1	2
			Spring, Arm Lock	1	2
<p>NOTE: A - Includes items 2-18. B - Includes items 17 and 18. C - Not shown. Includes items 6-9. D - Specify upholstery color when ordering. See page titled "Frame Finishes, Upholstery Colors, Decals, Literature and Packaging".</p>					