



9000XT, 9000XDT - Removable Space Saver Fixed Height Arms



9000XT, 9000XDT - Removable Space Saver Fixed Height Arms

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	8881130848	Kit, Armrest, Desk Length, Color - Right (Style 88)	1	1
2			Tube, Armrest, Spacer Saver, Desk, Color - Right	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Desk Length - Right	1	1
9			Arm Pad, Desk Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	1
11			Screw, Phillips Pan Head Tap (#6-20 x 3/8")	1	2
3		1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7	A,B	1115269	Guard, Clothing, Desk Length - Right	1	1
1		8881130849	Kit, Armrest, Desk Length, Color - Left (Style 88)	1	1
2			Tube, Armrest, Spacer Saver, Desk, Color - Left	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Desk Length - Left	1	1
9			Arm Pad, Desk Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
11			Screw, Phillips Pan Head Tap (#6-20 x 3/8")	1	2
3	B,E	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115270	Guard, Clothing, Desk Length - Left	1	1
1		8881130846	Kit, Armrest, Desk Length, Black - Right (Style 88)	1	1
2			Tube, Armrest, Spacer Saver, Desk, Black - Right	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Desk Length - Right	1	1
9			Arm Pad, Desk Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
11			Screw, Phillips Pan Head Tap (#6-20 x 3/8")	1	2
3	B,E	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115269	Guard, Clothing, Desk Length - Right	1	1
1		8881130847	Kit, Armrest, Desk Length, Black - Left (Style 88)	1	1
2			Tube, Armrest, Spacer Saver, Desk, Black - Left	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Desk Length - Left	1	1
9			Arm Pad, Desk Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
11			Screw, Phillips Pan Head Tap (#6-20 x 3/8")	1	2

9000XT, 9000XDT - Removable Space Saver Fixed Height Arms

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
3	B,E	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115270	Guard, Clothing, Desk Length - Left	1	1
1		8881130850	Kit, Armrest, Full Length, Black - Right (Style 78)	1	1
2			Tube, Armrest, Space Saver, Full Length, Black - Right	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Full Length - Right	1	1
9			NLA - Pad, Arm, Full Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
3	B,E	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115271	Guard, Clothing, Full Length - Right	1	1
1		8881130851	Kit, Armrest, Full Length, Black - Left (Style 78)	1	1
2			Tube, Armrest, Space Saver, Full Length, Black - Left	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Full Length - Left	1	1
9			NLA - Pad, Arm, Full Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
3	A,B	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115272	Guard, Clothing, Full Length - Left	1	1
1		8881130852	Kit, Armrest, Full Length, Color - Right (Style 78)	1	1
2			Tube, Armrest, Space Saver, Full Length, Color - Right	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Full Length - Right	1	1
9			NLA - Pad, Arm, Full Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
3	A,B	1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115271	Guard, Clothing, Full Length - Right	1	1
1		8881130853	Kit, Armrest, Full Length, Color - Left (Style 78)	1	1
2			Tube, Armrest, Space Saver, Full Length, Color - Left	1	1
3			Screw, Self Tapping (#6-20 x 1/2")	1	4
4			Handle, Quad Release	1	2
5			Spring, Torsion with Button Lock	1	2
6			Plug, Armrest, Black (3/4") (18GA)	1	2
7			Guard, Clothing, Full Length - Left	1	1
7			Guard, Clothing, Full Length - Left	1	1

9000XT, 9000XDT - Removable Space Saver Fixed Height Arms

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
9			NLA - Pad, Arm, Full Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
3		1116900	Screw, Self Tapping (#6-20 x 1/2")	1	4
7		1115272	Guard, Clothing, Full Length - Left	1	1
8	C,E	8881127509	Kit, Arm Pad w/Screws - Full Length	1	1
9			NLA - Pad, Arm, Full Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
8	C,E	8881127510	Kit, Arm Pad w/Screws - Desk Length	1	1
9			Arm Pad, Desk Length w/ Law Label	1	1
10			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	2
12	D	1058139	Kit, Armrest Hardware	1	1
13			NLA - Screw, Hex Head (1/4-20 x 1-3/8")	1	1
14			Cup, Socket Pivot	1	1
15			Locknut (1/4-20)	1	1
NOTE: A - Specify frame finish and upholstery color. See page titled "Frame Finishes, Upholstery Colors, Decals and Packaging." B - Includes items 2-7, 9 and 10 for Full Length Arms and items 2-7 and 9-11 for Desk Length Arms C - Includes items 9 and 10. D - Includes items 13-15. E - Specify upholstery color. See page titled "Frame Finishes, Upholstery Colors, Decals and Packaging."					