


MK5 SPJ-80 Upgrade to MK6

For wheelchairs with MK5 electronics and power tilt that need upgraded to MK6i electronics, the wheelchair will need to be returned to Invacare Tech Repair for the upgrade. Call Invacare Customer Service at 1-800-333-6900 for assistance.



MK5 SPJ-80 Upgrade to MK6

For wheelchairs with MK5 electronics and power tilt that need upgraded to MK6i electronics, the wheelchair will need to be returned to Invacare Tech Repair for the upgrade. Call Invacare Customer Service at 1-800-333-6900 for assistance.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1164850	Kit, MK5 SPJ-80 Upgrade to MK6 w/o Tilt	1	1
2			Joystick, SPJ+	1	1
3			NLA - Controller, MK5 NX-75	1	1
4			Cable, Joystick Extension	1	1
5			Cable, Joystick SPJ+	1	1
6			Locknut (#8-32)	1	2
7			Screw, Phillips Pan Head (#8-32 x 1")	1	2
8			Tie Wrap (11-1/2" L)	1	4
			NLA - Programming Guide, MK6i	1	1
			Instruction Sheet, MK5 SPJ-80 Upgrade to MK6	1	1

NOTE: A - Includes items 2-8, programming guide and instruction sheet.