


Hublock Replacement Kits

Attendant Foot Operated



Hublock Replacement Kits Attendant Foot Operated

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1164343	Kit, Hublock Cables and Hardware, Foot Operated	1	1
2			Locknut, Thin (1/4-20)	1	6
3			NLA - Washer, Steel (1/4 x 1/2 x 3/64")	1	4
4			Latch with Cable Assembly, Hublock - Left	1	1
4			Latch with Cable Assembly, Hublock - Right	1	1
5			NLA - Bracket, Cable/Wire Mounting, Black	1	2
6			Actuator Assembly, Foot Operated	1	1
7			Screw, Socket Button Head (1/4-28 x 1-1/4")	1	2
8			Tie Wrap (5-5/8" L)	1	3
4	B	1132358	Latch with Cable Assembly, Hublock - Left	1	1
4	B	1132357	Latch with Cable Assembly, Hublock - Right	1	1

NOTE: A - Includes items 2-8
 B - it is recommended to order both Hublock Latch with Cable Assembly part numbers 1132357 and 1132358 together to allow even wear on the components.