




# Casters and Mounting Hardware




Caster Models 3SOFR,  
CRFA654, CRFA655




Caster Model CRFF3




Caster Model CRF884



Caster Model CRF654



Caster Model CRF655



## Casters and Mounting Hardware

| Item Number                                                                                                                                                                                | Note | Part Number | Description                                            | Quantity Per |          |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|-------------|--------------------------------------------------------|--------------|----------|
|                                                                                                                                                                                            |      |             |                                                        | Package      | Assembly |
|                                                                                                                                                                                            |      | See Chart   | Caster Models                                          | 1            | 1        |
| 1                                                                                                                                                                                          |      | 1157282-NLA | NLA - Wheel Assembly, Soft Roll (3" x 1-1/2")          | 1            | 1        |
| 1                                                                                                                                                                                          |      | 1128070     | Wheel Assembly, Soft Roll (4" x 1")                    | 1            | 1        |
| 1                                                                                                                                                                                          |      | 1128071     | Wheel Assembly, Soft Roll (5" x 1-1/2")                | 1            | 1        |
| 2                                                                                                                                                                                          |      | 1032190-NLA | NLA - Wheel Assembly (3" x 1")                         | 1            | 1        |
| 3                                                                                                                                                                                          |      | 1060910     | Wheel Assembly, Composite Urethane, Light Up (4" x 1") | 1            | 1        |
| 4                                                                                                                                                                                          |      | 1114142     | Wheel Assembly, Semi-Pneumatic, Black (4" x 1")        | 1            | 1        |
| 5                                                                                                                                                                                          | A    | 1164311     | Wheel Assembly, Composite, Urethane (5" x 1")          | 1            | 1        |
| 6                                                                                                                                                                                          |      | 1159772     | Kit, Front Wheel Hardware                              | 1            | 1        |
|                                                                                                                                                                                            |      |             | Spacer (5/16 x 7/16)                                   | 1            | 2        |
|                                                                                                                                                                                            |      |             | Locknut (5/16-24)                                      | 1            | 1        |
|                                                                                                                                                                                            |      |             | Screw, Hex Head (5/16-24 x 2-1/4")                     | 1            | 1        |
| NOTE: A - After 10/4/09, the 5SOFT (5") Semi-Pneumatic is no longer available and is replaced by the 5SOFT (5") Urethane (1164311). Order in pairs for chairs manufactured before 10/5/09. |      |             |                                                        |              |          |

| Caster Model | Type                       | Part Number |
|--------------|----------------------------|-------------|
| 3SOFR        | Aluminum Soft Roll Casters | 1157282-NLA |
| CRFA654      | Aluminum Soft Roll Casters | 1128070     |
| CRFA655      | Aluminum Soft Roll Casters | 1128071     |
| CRFF3        | Micro Casters              | 1032190-NLA |
| CRF884       | Lighted Front Wheel        | 1060910     |
| CRF654       | Urethane Casters           | 1114142     |
| CRF655       | Urethane Casters           | 1164311     |