

Quick Reference Guide

Bed Model and Component Numbers

Manual

Split Spring Models

Fixed Height	5107-NLA	5880M004	5000-NLA	5770-NLA	7001M012	5100-NLA	5300M053-N LA
Two Crank Hi/Lo	5207-NLA	5880M004	5000-NLA	5770-NLA	7001M012	5200-NLA	5200M026
Single Crank Hi/Lo -U.S. Models	5307-NLA	5880M004	5000-NLA	5770-NLA	7001M012	5301-NLA	5300M053-N LA
Single Crank Hi/Lo -International Models	5302	5880M004	5000-NLA	5772	5770M014-N LA	5301-NLA	5300M053-N LA
Single Crank Hi/Lo	5307IVC	5880M004	5000IVC	5770IVC	7001M012	5301IVC	1116636-NLA See Note(s) A
Single Crank Hi/Lo	5307IVC	5880M004	5000IVC	5770IVC	7001M012	5301IVC	1157942 See Note(s) A
Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Cranks	Bed Ends	Hi/Lo Crank

NOTE: A - All models designated for "Government" use 1157942. All others use 1116636.

Semi-Electric

Split Spring Models

Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Cranks	Bed Ends	Hi/Lo Crank
Non-Synchro nizer (Before 7/90) Fixed Height	5108-NLA	5880M004	5000-NLA	5887-NLA	7001M012	5100-NLA	5300M053-N LA
Non-Synchro nizer (Before 7/90) Two Crank Hi/Lo	5208-NLA	5880M004	5000-NLA	5887-NLA	7001M012	5200-NLA	5200M026
Non-Synchro nizer (Before 7/90) Single Crank Hi/Lo	5308-NLA	5880M004	5000-NLA	5887-NLA	7001M012	5300-NLA	5300M053-N LA
Synchronizer (After 6/90) Fixed Height	5110-NLA	5880M004	5000-NLA	5890	7001M012	5100-NLA	5300M053-N LA
Synchronizer (After 6/90) Two Crank Hi/Lo	5210-NLA	5880M004	5000-NLA	5890	7001M012	5200-NLA	5200M026

Quick Reference Guide Bed Model and Component Numbers

Semi-Electric

Split Spring Models

Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Cranks	Bed Ends	Hi/Lo Crank
Synchronizer (After 6/90) Single Crank Hi/Lo - US Models	5310-NLA	5880M004	5000-NLA	5890	7001M012	5301-NLA	5300M053-NLA
Synchronizer (IVC Models) Single Crank Hi/Lo - US Models	5310IVC	5880M004	5000IVC	5890IVC	N/A	5301IVC	1116636-NLA See Note(s) A
Synchronizer (IVC Models) Single Crank Hi/Lo - US Models	5310IVC	5880M004	5000IVC	5890IVC	N/A	5301IVC	1157942 See Note(s) A
Synchronizer (After 6/90) Single Crank Hi/Lo - Canadian Models	5311-NLA	5880M004	5000-NLA	5891	7001M012	5301-NLA	5300M053-NLA
Synchronizer (After 6/90) Single Crank Hi/Lo - International Models	5312	5880M004	5000-NLA	5892-NLA	7001M012	5301-NLA	5300M053-NLA
Synchronizer (IVC Models) Single Crank Hi/Lo - Canadian Models	5311IVC	5880M004	5000IVC	5891IVC	N/A	5301IVC	1116636-NLA See Note(s) A
Synchronizer (IVC Models) Single Crank Hi/Lo - Canadian Models	5311IVC	5880M004	5000IVC	5891IVC	N/A	5301IVC	1157942 See Note(s) A

NOTE: A - All models designated for "Government" use 1157942. All others use 1116636.

Quick Reference Guide Bed Model and Component Numbers

Full Electric

Split Spring Models

Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Crank	Bed Ends	Hi/Lo Crank
Non-Synchronizer (Before 7/90)	5409-NLA	5880M004	5000-NLA	5997-NLA	5770M014-N LA	5400-NLA	1134790-NLA See Note(s) A
Synchronizer - U.S. Models (After 6/90) Homecare Beds	5410-NLA	5880M004	5000-NLA	5490-NLA	5770M014-N LA	5401-NLA	1134790-NLA See Note(s) A
Synchronizer - U.S. Models (IVC Models) Homecare Beds	5410IVC	5880M004	5000IVC	5490IVC	5770M014-N LA	5301IVC	1116636-NLA See Note(s) B
Synchronizer - U.S. Models (IVC Models) Homecare Beds	5410IVC	5880M004	5000IVC	5490IVC	5770M014-N LA	5301IVC	1157942 See Note(s) B
Synchronizer - U.S. Models (After 6/90) Bariatric Bed	BAR600	5880M004	BAR5000	BAR5490	5770M014-N LA	BAR5401	1134790-NLA See Note(s) A
Synchronizer - U.S. Models (IVC Models) Bariatric Bed	BAR600IVC	5880M004	BAR5000IVC	BAR5490IVC	N/A	BAR5301IVC	N/A
Synchronizer - Canadian Models	5411-NLA	5880M004	5000-NLA	5491	5770M014-N LA	5401-NLA	1134790-NLA
Synchronizer (IVC Models) - Canadian Models	5411IVC	5880M004	5000IVC	5491IVC	5770M014-N LA	5301IVC	1116636-NLA See Note(s) B
Synchronizer (IVC Models) - Canadian Models	5411IVC	5880M004	5000IVC	5491IVC	5770M014	5301IVC	1157942 See Note(s) B
Synchronizer - International Models	5412	5880M004	5000-NLA	5492-NLA	5770M014-N LA	5401-NLA	1134790-NLA See Note(s) A

Quick Reference Guide Bed Model and Component Numbers

Full Electric

Split Spring Models

Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Crank	Bed Ends	Hi/Lo Crank
Hi/Lo Crank Direct Replacement for Non-IVC style beds	N/A	N/A	N/A	N/A	N/A	N/A	1134790-NLA See Note(s) A
NOTE: A - Includes the Gearbox and the 80" Drive Shaft. B - All models designated for "Government" use 1116636. All others use 1157942.							

Quick Reference Guide Bed Model and Component Numbers

Careguard by Invacare Single Drive

Split Spring Models

Bed Style	Model Number	Emergency Handle	Head Springs	Foot Springs	Head/Foot Crank	Bed Ends	Hi/Lo Crank
Single Drive Crank Hi/Lo	CGBED-NLA	5880M004	5000-NLA	8050DK-NLA	5770M014-N LA	5301-NLA	1134790 See Note(s) A

NOTE: A - Includes the Gearbox and the 80" Drive Shaft.