






# Wheel Locks & Attaching Hardware



## Wheel Locks & Attaching Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1130628	Wheel Lock Assembly, Pull to Lock - Right	1	1
1	A	1130629	Wheel Lock Assembly, Pull to Lock - Left	1	1
2	A	1130630	Wheel Lock Assembly, Push to Lock - Right	1	1
2	A	1130631	Wheel Lock Assembly, Push to Lock - Left	1	1
3	B	1097058-NLA	NLA - Hill Holder Wheel Lock Assembly - Right	1	1
3	B	1097059-NLA	NLA - Hill Holder Wheel Lock Assembly - Left	1	1
4	C	1068210-NLA	Kit, Wheel Lock Attaching Hardware (1" Tubing)	1	1
5			Screw, Allen Head (1/4-20 x 1")	1	2
6			Clamp, Half, Non-Threaded	1	1
7			Clamp, Half, Threaded	1	1
8	B	1093376-NLA	NLA - Wheel Lock Assembly, Undermount - Right	1	1
8	B	1093377-NLA	NLA - Wheel Lock Assembly, Undermount - Left	1	1
9	D	1160548	Clamp Assembly, Wheel Lock	1	1
10		1130584	Package, Cap, Brake Lever End, Black	10	1
11		1075507	Tip, Plastic, Black	1	1
12	E	1146498	Handle, Wheel Lock Extension, Black (6") (Pair)	1	1

NOTE: A - 1 each of Item 10 included with wheel lock  
 B - 1 each of Item 11 included with wheel lock  
 C - For use with push, pull and hill holder wheel lock.  
 D - For use with undermount wheel lock.  
 E - Wheel Lock Extension Handle used with item number 1 and 2 only