


Mattress Assembly - 3500S



Mattress Assembly - 3500S

Standard Mattress Width

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9170101-NLA	NLA - Pad, Therapy	1	1
2	B	9170021-NLA	NLA - Sack, Air (28")	1	3
3		9170033-NLA	NLA - Mat, Rubber (82" x 36")	1	1
4		9172032-NLA	NLA - Retainer, Hook and Loop	1	1
5		9540040	NLA - Hose, Foot Section, Green	1	1
6		9540038	NLA - Hose, Seat Section, Red	1	1
7		9540036	NLA - Hose, Back Section, Yellow	1	1
8		9172027-NLA	NLA - Flag with Hardware, CPR	1	1
	A	BB9572038	NLA - Mattress Assembly, Standard - 3500S	1	1

NOTE: A - Not Shown
B - Item 4 must be ordered when replacing this part; it cannot be taken off with the old and used with the new.

Narrow Mattress Width

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		9170104-NLA	NLA - Pad, Therapy	1	1
2		9170093-NLA	NLA - Sack, Air (28")	1	3
3		9170070-NLA	NLA - Mat, Rubber (82" x 32")	1	1
4		9172032-NLA	NLA - Retainer, Hook and Loop	1	1
5		9540040	NLA - Hose, Foot Section, Green	1	1
6		9540038	NLA - Hose, Seat Section, Red	1	1
7		9540036	NLA - Hose, Back Section, Yellow	1	1
8		9172027-NLA	NLA - Flag with Hardware, CPR	1	1
	A	BB9572040-NLA	NLA - Narrow 3500S Mattress Assembly - Complete	1	1

NOTE: A - Not Shown