

Rear Wheels

Wheel Style	Single Wheel Assembly	Tire	Tube/Insert	Bearing	Bearing Spacer	Plugs
Composite Pneumatic (20 x 1-3/8") 1763BK	1034126 See Note(s) A, B	1017519	1017520	40117X024	1015756	N/A
Composite Pneumatic w/ Flat Free Insert (20 x 1-3/8") 4763BK	1034127 See Note(s) A, B	1017519	1023322	40117X024	1015756	1030936 See Note(s) C
Composite w/ Urethane (20 x 1-1/4") U1763BK	1038555 See Note(s) B, A	1037774	N/A	40117X024	1015756	N/A
Composite Pneumatic (22 x 1-3/8") 1764BK	1027149 See Note(s) A, B	1018544-NLA	1018545	40117X024	1015756	N/A
Composite Pneumatic w/ Flat Free Insert (22 x 1-3/8") 4764BK	1027150 See Note(s) A, B	1018544-NLA	1023323-NLA	40117X024	1015756	1030936 See Note(s) C
Composite w/ Urethane (22 x 1-1/4") U1764BK	1032115 See Note(s) A, B	1028250	N/A	40117X024	1015756	N/A
Composite DTM (24 x 7/8")	1080378 See Note(s) A, D	1080400 See Note(s) D	N/A	40117X024	1019695-NLA	N/A
Composite Corded Urethane (24 x 1") U2222C	1122726 See Note(s) A	1118209 See Note(s) E	N/A	40117X024	1019695-NLA	N/A
Composite Pneumatic (24 x 1-3/8") 2244BK	1025123 See Note(s) A	40176X060	40176X061-NL A	40117X024	1015756	N/A
Composite Pneumatic w/ Flat Free Insert (24 x 1-3/8") 4244BK See Note(s) C	1025124 See Note(s) A	40176X060	1010906	40117X024	1015756	1030936 See Note(s) C

Rear Wheels

[illegible]