

## Rear Wheels

Wheel Style	Composite 1763BK	Composite 4763BK	Composite U1763BK	Composite 2244BK	Composite 2222BK	Composite	Composite 4244BK
Wheel Size	20"	20"	20"	24"	24"	24"	24"
Single Tire Type	Pneumatic (20 x 1-3/8")	Pneumatic w/Flat Free Insert (20 x 1-3/8")	Urethane (20 x 1-1/4")	Urethane (24 x 1-1/4")	DTM (24 x 7/8")	Pneumatic (24 x 1-3/8")	Pneumatic w/Flat Free Insert (24 x 1-3/8")
Wheel Assembly	1034126 See Note(s) A	1034127 See Note(s) A	1038555 See Note(s) A	1032117 See Note(s) A	1080378 See Note(s) A, B	1025123 See Note(s) A	1025124 See Note(s) A
Wheel Assembly	N/A	N/A	N/A	N/A	1122726 See Note(s) A, D	N/A	N/A
Tire	1017519	1017519	1037774	1026904	1080400 See Note(s) B	40176X060	40176X060
Tire	N/A	N/A	N/A	N/A	1118209 See Note(s) C	N/A	N/A
Tube/Insert	1017520	1023322	N/A	N/A	N/A	40176X061	1010906
Hub Bearing	40117X024	40117X024	40117X024	40117X024	40117X024	40117X024	40117X024
Spacer	1015756	1015756	1015756	1015756	1019695-NLA	1015756	1015756
<p>NOTE: A - Includes tire, tube/insert, bearing and spacer. Does not include handrim.            B - Non-Corded Tire replaced with a Corded Tire as of 11/16/09. Refer to part #1122726 for a wheel assembly or part #1118209 for a tire. It is recommended to order (2) each in quantity to match.            C - Corded Tires maybe more difficult to stretch over the rim. The complete wheel assembly that has the tire and wheel rim preassembled maybe advisable.            D - The wheel assembly is now a U2222C Composite Corded Urethane 24 x1"</p>							