


Crutches

Forearm - Model 6152, 6153, 6154

Model 6152 - Junior

Model 6153 - Adult

Model 6154 - Tall Adult



Single Side Shown

Crutches

Forearm - Model 6152, 6153, 6154

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1049537-NLA	NLA - Cuff and Support Sub-Assembly, Double Snap, Black	1	1
1	A	1095678-NLA	NLA - Cuff and Support Sub-Assembly, Single Snap, Grey	1	1
1	B	1088819-NLA	NLA - Cuff and Support Sub-Assembly, Double Snap, Grey	1	1
2		1049536-NLA	NLA - Button, Snap, Double and Spring Wire	1	1
3		1049538-NLA	NLA - Tip, Rubber, Black (7/8" Tubing)	1	1
3		6135-NLA	NLA - Kit, Tips Extra Large Rubber, Grey (2" O.D. Base) Tip, Crutch, Grey - Extra Large	1 2	1 1
NOTE: A - For models built before 8/1/99. B - For models built after 8/1/99.					