


Control Panel - Front View



Control Panel - Front View

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
2		2000507	Palnut (9/16-18)	1	3
3		1110824-NLA	NLA - Tech Only - Circuit Breaker (4 AMP)	1	1
4		2002939-NLA	NLA - Tech Only - Hour Meter Assembly with Clip	1	1
5		1139950	Switch, Powe On/Off	1	1
6		1121475-NLA	NLA - Label, Control Panel, English - IRC5LX02AW	1	1
6		1121476	Label, Control Panel, Spanish - IRC5LX02AW	1	1
6		1121477	Label, Control Panel, French - IRC5LX02AW	1	1
6		1121478	Label, Control Panel, German - IRC5LX02AW	1	1
6		1121479	Label, Control Panel, Port/Brazil - IRC5LX02AW	1	1
6		1121480	Label, Control Panel, Norwegian - IRC5LX02AW	1	1
6		1121481	Label, Control Panel, Finland - IRC5LX02AW	1	1
6		1121482	Label, Control Panel, Danish - IRC5LX02AW	1	1
6		1142731	Label, Control Panel, Italian - IRC5LX02AW	1	1
6		1142732	Label, Control Panel, Polish - IRC5LX02AW	1	1
6		1142733	Label, Control Panel Text, Czech - IRC5LX02AW	1	1
6		1142734	Label, Control Panel, Greek - IRC5LX02AW	1	1
6		1142735	Label, Control Panel, Dutch - IRC5LX02AW	1	1
6		1142736	Label, Control Panel, Swedish - IRC5LX02AW	1	1
6		1164897	Label, Control Panel, Hungarian - IRC5LX02AW	1	1
7	C	1110561	Flowmeter without Palnuts - 5 Liter	1	1
8		2000473	Connector, Barbed (1/4")	1	1
9		1140334-NLA	NLA - Decal, Logo - Platinum S AWQ	1	1
10	A	1128379-NLA	NLA - Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
10	B	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
11		1101591	Support, Control Panel	1	2
12		1112562-NLA	NLA - Tech Only - Seal, Circuit Breaker	1	1

NOTE: A - Knob is Silver with Black increasing flow symbol. Use on flowmeters with Phillips Head Indented Screw on top.
 B - Knob is Black with Silver increasing flow symbol. Use on flowmeters with Slotted Head Screw on top.
 C - Palnuts are not included, refer to Item #2.