



Heel Loops



Heel Loop



Heel Loop with Ankle Strap

For Adjustable Angle Footplates - AT5543

Footplate Size	Chair Width	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Single # Right	Heel Loop w/ Ankle Strap Single # Left	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
Medium (4-1/4" Wide)	14-15" Wide or Tapered 16" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1035931
Large (5-1/4" Wide)	16-17" Wide or Tapered 17-18" Wide	1017621	1035917-NLA	1053313-NLA	1053314-NLA	1053329	1035931
X-Large (6-1/4" Wide)	18-22" Wide or Tapered 19-22" Wide	1027000	1035929	1035932-NLA	1035933-NLA	1035930	1035931
XX-Large (7-1/4" Wide)	18-22" Wide or Tapered 19-22" Wide	1027000	1035929	1035932-NLA	1035933-NLA	1035930	1035931

Heel Loops

For Composite Footplates - 1651

Footplate Size	Chair Width	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Single # Right	Heel Loop w/ Ankle Strap Single # Left	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
Small (6" Wide)	14" Wide or Tapered 16-17" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1035931
Medium (7" Wide)	15-16" Wide or Tapered 18-22" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1035931
Large (8" Wide)	17-22" Wide or Tapered 19-22" Wide	1027000	1035929	1035932-NLA	1035933-NLA	1035930	1035931

Heel Loop Attaching Hardware Kit

1035931
See Note(s) A
(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")
NOTE: A - Includes hardware for (1) pair of footplates