


Walker (Single Release)

Model 6281-A (4 Pack)

Model 6281-I (Adult Single Pack)

Model 6281-JR (4 Pack)



Walker (Single Release)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	1074462-NLA	NLA - Button, Snap	1	4
1	D	1130780-NLA	NLA - Kit, Snap Buttons	10	1
1		1094945	Button, Snap	1	4
2		1040655	Anti-Rattle, Grey (1-1/8")	1	4
3	A,B	1058000	Leg Extension Assembly (10")	1	2
4		1040656	Tip, Rubber, Dark Grey (1-1/8")	1	4
5		1060558-NLA	NLA - Folding Mechanism - Adult	1	1
5		1060559-NLA	NLA - Folding Mechanism - Junior	1	1
6		17564X000-NLA	NLA - Screw, Phillips Pan Head (#10-24 x 1-3/4")	1	4
<p>NOTE: A - Due to style changes, when replacing Leg Extensions before 6/25/98, order (4) each. Models built after 6/24/98 can order Single Assemblies. B - Includes Leg and Items 2 and 4. C - Individual item no longer available, use 1130780 as its replacement. D - Kit is No Longer Available, use 1094945 as the replacement.</p>					