

# Rollator (Four Wheeled) - Model 65851R - 65851B

Available in Red (65851R) and Blue (65851B)



# Rollator (Four Wheeled) - Model 65851R - 65851B

Available in Red (65851R) and Blue (65851B)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1151023	Grip, Hand (Pair)	2	1
2		1165495-NLA	NLA - Brake Assembly w/ Handle and Grip - Right	1	1
3		1165496-NLA	NLA - Brake Assembly w/ Handle and Grip - Left	1	1
4		1151020-NLA	NLA - Backrest Assembly, Red	1	1
4		1165096-NLA	NLA - Backrest Assembly, Blue	1	1
5	B,F	1151021-NLA	NLA - Knob w/ Hardware (Pair)	2	1
5	G	1169007	Kit, Adjustable Knob (Pair)	2	1
6		1151015-NLA	NLA - Seat w/ Mounting Hardware	1	1
7		1151014-NLA	NLA - Basket	1	1
8		1151024-NLA	NLA - Package, Cap	6	1
9	C	1151022-NLA	NLA - Package, Springs, Brake (Pair)	2	1
10	D	1151019-NLA	NLA - Kit, Caster with Hardware - Rear (Pair)	1	1
11	E	1151018-NLA	NLA - Kit, Caster with Fork and Hardware - Front (Pair)	1	1

NOTE: A - Includes (2) Hand Grips.  
 B - Includes (2) Knobs and Attaching Hardware.  
 C - Includes (2) Springs.  
 D - Includes (2) Rear Casters with Attaching Hardware.  
 E - Includes (2) Front Casters and Fork with Attaching Hardware.  
 F - Use with Rollators manufactured before 02/06/12.  
 G - Use with Rollators manufactured after 02/05/12.

Specify Model Color when ordering.