



Walker (Dual Release)

NOTE: 6291 Series Walkers Side Cross Brace was changed from plastic to aluminum after 8/19/17.

Model 6261 - A - Adult Walker (4 Pack)
Model 6261 - JR (4 Pack)



Model 6291 - A (4 Pack)
Model 6291 - I (Adult Single Pack)
Model 6291 - JR (4 Pack)



Walker (Dual Release)

NOTE: 6291 Series Walkers Side Cross Brace was changed from plastic to aluminum after 8/19/17.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	D	1074462-NLA	NLA - Button, Snap	1	4
1	F	1130780-NLA	NLA - Kit, Snap Buttons	10	1
1		1094945	Button, Snap	1	4
2		1040655	Anti-Rattle, Grey (1-1/8")	1	4
3		1040656	Tip, Rubber, Dark Grey (1-1/8")	1	4
4	A,B	1058000	Leg Extension Assembly (10")	1	4
5	C,E	1043538-NLA	NLA - Nut, Push	1	4
5		1130784-NLA	NLA - Kit, Push Nut	10	1

NOTE: A - Due to style changes, when replacing Leg Extensions before 6/25/98, Order (4) each. Models built after 6/24/98 can order Single Assemblies.
 B - Includes Leg and Items 2 and 3.
 C - To be used with Models 6291-A, 6291-1 and 6291 Jr.
 D - Individual item no longer available, use 1130780 as its replacement.
 E - Individual item no longer available, use 1130784 as its replacement.
 F - Kit is no longer available, use 1094945 as its replacement.