

# **A: READ THIS FIRST (THIS CATALOG IS DIFFERENT!)**

TDX SI and TDX SI-HD  
parts catalog information is now divided into two separate catalogs.

The Form Numbers and Cover Titles are listed below.

**Form # 08-080**  
TDX SI and TDX SI-HD  
Powered Wheelchair

**Form # 06-110**  
MK6i Electronics  
(MK6i was introduced on or about 10/9/06)

**Use both form numbers to order a “complete” Parts Catalog**

**NOTE:** TDX SI-2 model uses the MK5 NX-75 Controller and SPJ+ Joystick