


Seat Pans (After 4/18/10)



Seat Pans (After 4/18/10)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Seat Pans (12" - 16" Seat Depth Range)	1	1
1	A	See Chart	#1 Seat Pans (15" - 19" Seat Depth Range)	1	1
1	A	See Chart	#1 Seat Pans (18" - 22" Seat Depth Range)	1	1
	B	1125723	Kit, Hook Fastener w/ Adhesive Back	1	1

NOTE: A - When ordering a new seat pan, the Hook Fastener w/ Adhesive Back Kit must be ordered as well. The Seat Pans are steel and accommodates both flush and drop mount hardware
 B - Not Shown.

#1 Seat Pans (12" - 16" Seat Depth Range)

Seat Depth Range	12" - 13" Wide	14" - 15" Wide	16" - 17" Wide	18" Wide
12" - 16" Deep	1165788	1165789	1165790	1165791

#1 Seat Pans (15" - 19" Seat Depth Range)

Seat Depth Range	12" - 13" Wide	14" - 15" Wide	16" - 17" Wide	18" - 19" Wide	20" - 21" Wide	22" - 23" Wide	24" Wide
15" - 19" Deep	1165792	1165793	1165794	1165795	1165796	1165797	1165798

#1 Seat Pans (18" - 22" Seat Depth Range)

Seat Depth Range	12" - 13" Wide	14" - 15" Wide	16" - 17" Wide	18" - 19" Wide	20" - 21" Wide	22" - 23" Wide	24" Wide
18" - 22" Deep	N/A	1165799	1165800	1165801	1165802	1165803	1165804