



Adjustable Height Arms

Flip Back FADJDESK and FADJFULL


Non-Flip Back ADJDESK and ADJFULL




Adjustable Arm Base - Non-Flip Back



Adjustable Arm Base - Flip Back



**HARMHT - High Height
Arm Base Assembly**



**SARMHT - Standard (Low) Height
Arm Base Assembly**

Adjustable Height Arms

Flip Back FADJDESK and FADJFULL

Non-Flip Back ADJDESK and ADJFULL

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881053578	Kit, Arm Pad w/ Screws, Desk Length (9" L x 2" W) (Pair)	2	1
1	A	8881053577	Kit, Arm Pad w/ Screws, Full Length (14" L x 2" W) (Pair)	2	1
2		1111714	Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	2	1
3		1165476	Arm Tube Assembly w/ Plug Buttons, Desk Length, Black - Right/Left	1	1
3		1165478	Arm Tube Assembly w/ Plug Buttons, Full Length, Black - Right/Left	1	1
4		1130582	Package, Plug Button, Black (7/8")	10	2
5	B	1159497	Arm Base Assembly, Black, ADJDESK/FULL - Right	1	1
			Base, Arm, Black - Right/Left	1	1
6			Spring, Arm Lock	1	1
7			Pin, Arm Lock, Black	1	1
8			Rivet	1	6
9			NLA - Guard, Clothing	1	1
10			Handle, Lock, Black - Right	1	1
11			Rivet, Tubular	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
5		1159498	Arm Base Assembly, Black, ADJDESK/FULL - Left	1	1
			Base, Arm, Black - Right/Left	1	1
6			Spring, Arm Lock	1	1
7			Pin, Arm Lock, Black	1	1
8			Rivet	1	6
9			NLA - Guard, Clothing	1	1
10			Handle, Lock, Black - Left	1	1
11			Rivet, Tubular	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
12	C	1159480	Arm Base Assembly, Black, FADJDESK, Full - Right (High)	1	1
13			Base, Arm, Black - Right/Left	1	1
15			Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
16			Rivet, Tubular	1	1
17			Handle, Lock, Black - Right	1	1
18			Spring, Arm Lock	1	1
19			Pin, Arm Lock, Black	1	1
20			Plug Button, Black (5/8")	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
12	C	1159481	Arm Base Assembly, Black, FADJDESK, Full - Left (High)	1	1
13			Base, Arm, Black - Right/Left	1	1
15			Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
16			Rivet, Tubular	1	1

Adjustable Height Arms

Flip Back FADJDESK and FADJFULL

Non-Flip Back ADJDESK and ADJFULL

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
17			Handle, Lock, Black - Left	1	1
18			Spring, Arm Lock	1	1
19			Pin, Arm Lock, Black	1	1
20			Plug Button, Black (5/8")	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
21	D	1136591	Arm Base Assembly, Black, FADJDESK, Full - Right (Standard)	1	1
14			Base, Arm, Black - Right/Left	1	1
15			Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
16			Rivet, Tubular	1	1
17			Handle, Lock, Black - Right	1	1
18			Spring, Arm Lock	1	1
19			Pin, Arm Lock, Black	1	1
20			Plug Button, Black (5/8")	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
21	D	1136592	Arm Base Assembly, Black, FADJDESK, Full - Left (Standard)	1	1
14			Base, Arm, Black - Right/Left	1	1
15			Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
16			Rivet, Tubular	1	1
17			Handle, Lock, Black - Left	1	1
18			Spring, Arm Lock	1	1
19			Pin, Arm Lock, Black	1	1
20			Plug Button, Black (5/8")	1	1
	E	1090768	Kit, Lock Handle w/Rivets (Pair)	1	1
NOTE: A - Specify upholstery color U67 when ordering B - Includes items 6 - 11 and Arm Base. C - Includes item 13 and items 15 - 20. Arm Base Assembly used for HARMHT High Arm height adjustments 12" - 16" in one inch increments. D - Includes items 14 thru 20. Arm Base Assembly used for SARMHT Standard (Low) height adjustments 9" - 13" in one inch increments. E - Includes both right and left Lock Handles and (2) ea. Tubular Rivets. Not shown in illustration.					