


# Shroud Covers, Recline Actuator Curve Contoura CGTR/CGTRE



## Shroud Covers, Recline Actuator Curve Contoura CGTR/CGTRE

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1166889	Shroud w/ Decal, Curve Contoura	1	1
3		1165887	Bracket, Recline Actuator Module, Black	1	1
4	A	1168999-NLA	NLA - Kit, Curve Contoura Recline Actuator Shroud Mounting Hardware	1	1
5			Bracket, Shroud Rear Mount Clevis, Black	1	1
6			Screw, Socket Head w/ Patch (5/16-18 x 2")	1	2
7			Washer, Flat (5/16 x 5/8 x 1/16")	1	2
8			Locknut (5/16-18)	1	2
9			Screw, Thumb (#10-32 x 1/2")	1	2

NOTE: A - Includes items 5-9

Item #2 Reflective Decal is not available for sale individually.