


# Back Post and Upholstery



## Back Post and Upholstery

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		440106-1	Back Post Assembly (15" W x 20" H)	--	1
1		440106-2	Back Post Assembly (16" W x 20" H)	--	1
1		440106-3	Back Post Assembly (17" W x 20" H)	--	1
1		440106-4	Back Post Assembly (18" W x 20" H)	--	1
1		440106-5	Back Post Assembly (19" W x 20" H)	--	1
1		440106-6	Back Post Assembly (20" W x 20" H)	--	1
1		1072266-1	Back Post Assembly (15" W x 24" H)	--	1
1		1072266-2	Back Post Assembly (16" W x 24" H)	--	1
1		1072266-3	Back Post Assembly (17" W x 24" H)	--	1
1		1072266-4	Back Post Assembly (18" W x 24" H)	--	1
1		1072266-5	Back Post Assembly (19" W x 24" H)	--	1
1		1072266-6	Back Post Assembly (20" W x 24" H)	--	1
2		440109	Post, Back - Right, AAB22 (22")	--	1
2		440110	Post, Back - Left, AAB22 (22")	--	1
2		440112	Post, Back - Right, AAB26 (26")	--	1
2		440113	Post, Back - Left AAB26 (26")	--	1
3	A	450009	Back Upholstery for widths (15" & 16")	--	1
3	A	450010	Back Upholstery for widths (17" & 18")	--	1
3	A	450011	Back Upholstery for widths (19" & 20")	--	1
4		1087493	Cable Assembly C45 (15 - 20" W)	--	1
5		1087595	Cap, Recessed	--	2
6		20228X013-NLA	NLA - Screw, Phillips Oval Head (#10 x 7/8")	--	2
7		R2763	Screw, Cross Pan Head (10-32 x 5/8")	--	1
8		1021260	Tie Wrap (5-5/8" L)	--	3
9		1032244	Grip, Foam Hand (9" L)	--	2
10		1043872-NLA	NLA - Bar, Spreader, Black (15" W)	--	1
		1043873	Bar, Spreader, Black (16" W)	--	1
		1043874	Bar, Spreader, Black (17" W)	--	1
		1043875	Bar, Spreader, Black (18" W)	--	1
		1043876	Bar, Spreader, Black (19" W)	--	1
		1043877	Bar, Spreader, Black (20" W)	--	1
11		450111	Hardware Back Mounting	--	1

NOTE: A - Specify width and height for item 3 and only available with old Style Back Post, item 1.