

# Miscellaneous Accessories

## Carrying Pocket

Width	Nylon	Vinyl
16 - 20" Wide	8881047959	8881047957 See Note(s) A

NOTE: A - Specify upholstery color when ordering. See page titled, "Frame Finishes, Touch-Up Paint and Upholstery Options."

## Auto Style Push Button Positioning Strap

Chair Width	Part Number	Seat/Chest Positioning Strap w/ Auto Style Buckle Length (Eyelet to Eyelet)
16" - 18" Wide	1004310	48" L
20" Wide & All Captain Seats	1034625	61" L
22" Wide	1034622-NLA See Note(s) A	69" L

NOTE: A - For replacement contact The Aftermarket Group (TAG).

## Airline Style Positioning Strap

Chair Width	Part Number	Seat/Chest Positioning Strap w/ Airline Style Buckle Length (Eyelet to Eyelet)
16" - 18" Wide	TAGTS1034751-NLA	46" L
20" Wide & All Captain Seats	1034752-NLA See Note(s) A	61" L
22" Wide Chairs	1034753-NLA	70" L


NOTE: A - 1034752 Airline Buckle Style is No Longer Available. Use Auto Style part number 1034625 as a replacement.

## Hook and Loop Style Positioning Strap


Chair Width	Part Number	Seat/Chest Positioning Strap w/ Hook and Loop Style Buckle Length (Eyelet to Eyelet)
16" - 18" Wide	1034744 See Note(s) A	45" L
20" Wide & All Captain Seats	1034745 See Note(s) A	61" L

NOTE: A - Hook and Loop Style Positioning Straps for power chairs used for chest positioning only.


## Miscellaneous Accessories



**Auto Style Push Button Buckle**



**Airline Style Buckle Positioning Strap**



**Hook and Loop Style Buckle Positioning Strap**