

H720T Glide About Commode



SPC_024A

H720T Glide About Commode

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1525864	Seat, Toilet w/ Holder, Gray	1	1
1		1525865	Seat, Toilet w/o Holder, Gray	1	1
2		1470947	Plate, Cover, Removable	1	1
2		1525866	Holder, Toilet Pan, Gray	1	1
3		1470948	Pan, Sanitary w/ Lid	1	1
4		1525868	Footrest, Right	1	1
4		1525867	Footrest, Left	1	1
5		1470952	Footplate, Left/Right	1	1
6		1525869	Armrest, Right, Gray	1	1
6		1525870	Armrest, Left, Gray	1	1
7		1470955	Pad, Armrest	1	2
8		1470956	Caster w/o Brake	1	2
9		1470957	Caster w/ Brake	1	2
10		1470958	Cover, Backrest	1	1
11		1470959	Knob, Backrest	1	2
12		1470960	Lock, Armrest	1	2
13		1470949	Tip, Rubber	1	2
14		1470962	Calf Strap	1	1