


# Rollator (Four Wheeled)

## Model 65500 - Stingray

Model 65500 (Stingray)  
Blue/Grape Finish



## Rollator (Four Wheeled) Model 65500 - Stingray

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1076180-NLA	NLA - Grip, Foam	2	1
2	A	1076184-NLA	NLA - Kit, Hand Brake	1	2
3	B	1076177-NLA	NLA - Kit, Back Support	1	1
4	C	1076182-NLA	NLA - Kit, Cane	1	2
5	D	1076175-NLA	NLA - Knob, Adjustable with Hardware	1	2
6		1076172-NLA	NLA - Seat	1	1
7	E	1076178-NLA	NLA - Seat Hardware	1	1
8	F	1076176-NLA	NLA - Kit, Wheel w/ Hardware	1	4
9		1076169-NLA	NLA - Basket	1	1
10		1076174-NLA	NLA - Tray	1	1
11		1082534-NLA	NLA - Cap, Plastic, Black	8	1
12		1094535-NLA	NLA - Fork, Front	1	2
13		1110750-NLA	NLA - Kit, Connecting Link w/ Hardware	1	2
14		1104754-NLA	NLA - Pad, Brake	1	2
<p>NOTE: A - Includes Brake Handle, Brake Holder and Brake Cable Only.            B - Includes Support Strap, Right and Left Anchor, Screw and Nuts.            C - Includes Grip, Hand Brake, Cable and Screws.            D - Includes Knob and Bolt.            E - Includes Brackets, Plates and Screws.            F - Includes Wheel, Bearings, Spacer, Nut, Screw and Caps.</p>					