



# Harnesses (continued)


As of 7/1/2009, switches are mercury free.




**Elevate Actuator and Manual Tilt Sensor Harness (Manual Tilt on TDX SPREE)**





**Elevate/Tilt Actuator and Sensor Harness (Formula CGT on TDX SPREE)**



**Tilt Actuator/MPCTMT Power Legrest and Sensor Harness (Formula CGT on TDX SC)**



**Tilt Actuator and Sensor Harness (Formula CGT, CGTE, CGTR, CGTRE on FDX)**



## Harnesses (continued)

As of 7/1/2009, switches are mercury free.

| Item Number   | Note | Part Number | Description  | Quantity Per |          |
|---|------|-------------|--|--------------|----------|
|   |      |             |  | Package      | Assembly |
| 1   | A    | 1151591-NLA | NLA - Harness, Elevate Actuator and Manual Tilt Sensor       | 1            | 1        |
| 2   | B    | 1151608-NLA | NLA - Harness, Elevate/Tilt Actuator and Sensor              | 1            | 1        |
| 3   | C    | 1153867-NLA | NLA - Harness, Tilt Actuator/MPCTMT Power Legrest and Sensor | 1            | 1        |
| 4   | D    | 1168442     | Harness, Tilt Actuator and Sensor                            | 1            | 1        |
| 5   |      | 1140099     | Switch, Drive Lockout Mercury Free Ball (MK6i) (50" L)       | 1            | 2        |
| 5   | E    | 1168753     | Switch, Drive Lockout Mercury Free Ball (MK6i) (18" L)       | 1            | 2        |
| 6   |      | 1122586     | Sensor, Proximity Switch (30 1/4" L)                         | 1            | 1        |
| 6   | F    | 1127290     | Sensor, Proximity Switch (7" L)                              | 1            | 1        |
| NOTE: A - Used for TDX SPREE with a Manual Tilt (MT) seating option<br>B - Used for TDX SPREE with a Formula CGT seating option<br>C - Used for TDX SC with a Formula CGT seating option with MPCTMT Power Legrest.<br>D - Used for FDX with a Formula CGT, CGTE, CGTR, CGTRE seating option<br>E - Used for FDX with a Formula CGTR/CGTRE seating option<br>F - Used for FDX with a Formula CGE seating option |      |             |  |              |          |