


I.V. Telescoping Rod

Not available for the Reclining Back (RECLINER) w/ Upholstery or the Recliner Back (RECLINEWU) w/o Upholstery



Base Frame, Seat Frame Assemblies
And Back Canes Are Shown For
Reference Only.

I.V. Telescoping Rod

Not available for the Reclining Back (RECLINER) w/ Upholstery or the Recliner Back (RECLINEWU) w/o Upholstery

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1159546	Kit, I.V., Telescoping Rod	1	1
2			I.V. Telescoping Rod Assembly, Chrome	1	1
3			Receiver, I.V. Telescoping Rod, Black	1	1
4			NLA - Screw, Hex Head (1/4-28 x 1-1/4")	1	1
5			NLA - Screw, Hex Head w/ Patch (1/4-28 x 3/4")	1	1
6			Screw, Hex Head (1/4-28 x 1-5/8")	1	2
7			Locknut (1/4-28)	1	3
8			Spacer, Coved (1/4 x 11/16 x 1/2")	1	2

NOTE: A - Includes items 2-8.