




Back Canes - Adjustable Height (ABH1720)



Back Canes - Adjustable Height (ABH1720)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1134189	Adjustable Height Cane Assembly (Pair)	1	1
2			Tube, Receiver Cane, Adjustable Height, Black	1	2
3			Tube, Adapter, Stroller Handle, Black	1	2
4			Screw, Allen Head (1/4-28 x 1-1/4")	1	6
5			Locknut (1/4-28)	1	6
6	B	1166079	Adjustable Height Cane Assembly (Pair)	1	1
8			Tube, Receiver Cane, Adjustable Height, Black	1	2
10			Handle, Push - Right	1	1
11			Handle, Push - Left	1	1
12			Grip, Hand, Foam (4")	1	2
13			Panel Fastener, Snap In	1	4
14			Screw, Socket Head (1/4-28 x 1-1/4")	1	4
15			Locknut, Thin (1/4-28)	1	4
7	C	1170184	Adjustable Height Cane Assembly (Pair)	1	1
9			Tube, Receiver Cane, Adjustable Height, Black	1	2
10			Handle, Push - Right	1	1
11			Handle, Push - Left	1	1
12			Grip, Hand, Foam (4")	1	2
13			Panel Fastener, Snap In	1	4
14			Screw, Socket Head (1/4-28 x 1-1/4")	1	4
15			Locknut, Thin (1/4-28)	1	4
16		1134187	Tube, Cane, Adjustable Height, Black (16 - 20")	1	2
NOTE: A - Includes items 2 - 5. Used only with Stroller Handle before 4/19/10. B - Includes items 8 and 10 - 15. Used with Standard Aluminum (SCANE), and Heavy Duty Package (HDPG) chairs. Also used with Rocker Back (DROCK) from 4/18/10 to 8/24/10. C - Includes items 9 - 15. Only used with Rocker Backs (DROCK).					