


Rollator (Four Wheeled - Discontinued 6/07)

Model 65700 - Frontier Lite

Model 65700 (Frontier Lite)
Burgundy Finish



Rollator (Four Wheeled - Discontinued 6/07) Model 65700 - Frontier Lite

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1076186	Grips, Hand	2	1
2		1099372-NLA	NLA - Cover	1	1
3		1094528-NLA	NLA - Pad, Backrest	1	1
4	A	1099632-NLA	NLA - Kit, Backrest Tube, Burgundy	1	1
5		1099367-NLA	NLA - Package, Anti-Rattle, Black (1")	2	1
6	B	1076175-NLA	NLA - Knob, Adjustable with Hardware	1	2
7		1099299-NLA	NLA - Seat	1	1
8	C	1099373-NLA	NLA - Kit, Seat Hardware	1	1
9	F	1110749-NLA	NLA - Kit, Connection Link w/ Hardware	1	2
10		1099369-NLA	NLA - Fork, Front	2	1
11	D	1099370-NLA	NLA - Kit, Wheel (5"), Front	2	1
12	E	1099371-NLA	NLA - Kit, Wheel - Rear (5")	2	1
13		1099368-NLA	NLA - Package, Rubber Tip	2	1
14		1098988-NLA	NLA - Basket	1	1
15		1076174-NLA	NLA - Tray	1	1
<p>NOTE: A - Includes Backrest Tube & Snap Buttons. B - Includes Knob and Bolt. C - Includes Brackets and Screws. D - Includes Wheel, Bearings, Bolt, Washer and Nut. E - Includes Wheel, Bearings, Bolt, Washer, Spacer and Nut. F - Includes Front and Rear Linkage, Hardware and Instruction Sheet.</p>					