


Sleep Surfaces - Link, Fabric



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		026636	Fabric, Link (76")	1	1
1		026637	Fabric, Link (80")	1	1
2		51000M136	Spring, Helical (18 turns)	1	12
3	A	51000M135	Spring, Helical (11 turns)	1	31
4	A	001477	Spring, Helical (11-3/4 Turn)	1	31
5		1113708-NLA	NLA - Bed Deck & Frame - 720 Style w/ 76" Link Fabric	1	1
5		1113709-NLA	NLA - Bed Deck & Frame - 720 Style w/ 80" Link Fabric	1	1
5		1113714-NLA	NLA - Bed Deck & Frame - SC900 Style w/ 76" Link Fabric	1	1
5		1113715-NLA	NLA - Bed Deck & Frame - SC900 Style w/ 80" Link Fabric	1	1
6		026510-NLA	NLA - Bed Frame - 720 Style w/ 76" Link Fabric	1	1
6		026513-NLA	NLA - Bed Frame - 720 Style w/ 80" Link Fabric	1	1
6		027051-NLA	NLA - Bed Frame - SC900 Style w/ 76" Link Fabric	1	1
6		027052-NLA	NLA - Bed Frame - SC900 Style w/ 80" Link Fabric	1	1
7		026518-NLA	NLA - Spring, Thigh Section	1	1
8		026521-NLA	NLA - Wire Deck, Foot Section	1	1
9		026947-NLA	NLA - Spring, Head (76")	1	1
9		026949-NLA	NLA - Spring, Head (80")	1	1
10		026159	Stop, Mattress	1	1

NOTE: A - 80" Link, Fabric requires an assembly quantity of 33 Helical Springs while the 76" Link, Fabric requires 31.