


Sleep Surfaces - Pan Deck



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	026763-NLA	NLA - Pan Deck Assembly, Head Section (76")	1	1
1	A	026764-NLA	NLA - Pan Deck Assembly, Head Section (80")	1	1
2		026503-NLA	NLA - Head Section (76")	1	1
2		026652-NLA	NLA - Head Section (80")	1	1
3		026506	Pan Deck Assembly (76")	1	1
3		026655-NLA	NLA - Assembly, Pan Deck (80")	1	1
4	B	026507	Pan Deck - Seat Section Assembly	1	1
5		026623	Rivet (1/4")	1	6
6	C	026766-NLA	NLA - Pan Deck Assembly, Thigh Section	1	1
7		026518-NLA	NLA - Spring, Thigh Section	1	1
8		026508	Pan Deck, Thigh Section	1	1
9	D	026767-NLA	NLA - Pan Deck Assembly - Foot Section	1	1
10		026521-NLA	NLA - Wire Deck, Foot Section	1	1
	E	026984-NLA	NLA - Kit, Sleep Surface Hardware	1	1
	E,F	000194	Includes: 4 ea. Nuts and 4 ea. Screws. Grommet, Pan Deck	1	6

NOTE: A - Includes: Items 2-3, 15 Rivets (Item #5) and Not Shown.
 B - Includes: Item 5.
 C - Includes: Items 7-8, 8 Rivets (Item #5) and Not Shown.
 D - Includes: Items 3, 10, 15 Rivets (Item #5) and Not Shown.
 E - Not Shown.
 F - 2 Grommets are used on each section (Head, Thigh, and Foot).