


# Sleep Surfaces - Wire Deck



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	026802-NLA	NLA - Wire Deck Assembly - Head Section (76")	1	1
1	A	026803-NLA	NLA - Wire Deck Assembly - Head Section (80")	1	1
2		026777-NLA	NLA - Wire Deck - Head Section (76")	1	1
2		026778-NLA	NLA - Wire Deck - Head Section (80")	1	1
3	B	026507	Pan Deck - Seat Section Assembly	1	1
4		026623	Rivet (1/4")	1	6
5	C	026805-NLA	NLA - Wire Deck Assembly - Thigh Section	1	1
6		026779-NLA	NLA - Wire Deck - Thigh Section	1	1
7	D	026806-NLA	NLA - Wire Deck Assembly - Foot Section	1	1
8		026521-NLA	NLA - Wire Deck, Foot Section	1	1
9		000046	Nut, Hex Head w/ Nylon Patch (1/4-20")	1	2
10		000124-NLA	NLA - Screw, Hex Head (1/4-20 x 3/4")	1	2
	E	026984-NLA	NLA - Kit, Sleep Surface Hardware Includes: 4 ea. Nuts and 4 ea. Screws.	1	1

NOTE: A - Includes: Item 2-4 and Not Shown.  
 B - Includes: Item 4.  
 C - Includes: Items 6 and Not Shown.  
 D - Includes: Items 8-10 and Not Shown.  
 E - Not Shown.