


# Bumpers


## Standard Carriage Discontinued 3/15/99



*Bed Frame shown  
for reference only*

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		012972-NLA	NLA - U-Bumper (76")	1	1
1		012971	U-Bumper (80")	1	1

## Maxi-Lock Bumper



*Bed Frame shown  
for reference only*

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	002886-NLA	NLA - U-Bumper w/ Screws (Dark Brown) (76")	1	1
1	A	1133068-NLA	NLA - U-Bumper w/ Screws (Taupe) (76")	1	1
1	A	002885-NLA	NLA - U-Bumper w/ Screws (Dark Brown) (80")	1	1
1	A	1125224-NLA	NLA - U-Bumper w/ Screws (Taupe) (80")	1	1
2		008682	Screw, Hex Head (1/4-20 x 5/8")	1	4

NOTE: A - These items are painted. Please make sure you select the correct color option.