


Bed Rails

Model IH5830 - 1/2 Length, Head End

Model IH5830 - Discontinued 2/05



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	027039-NLA	NLA - Side Rail Assembly (1/2 Length)	1	1
2		1113726-NLA	NLA - Side Rail, 1/2 Length - Right	1	1
2		1113727-NLA	NLA - Side Rail, 1/2 Length - Left	1	1
3		033786	Screw, Socket Head (#10-32 x 3/8")	1	4
	B	020258	Wrench, Allen (5/32")	1	1
	B	026715-NLA	NLA - Instruction Sheet, Bed Rail	1	1
	B	026197-NLA	NLA - Hardware, Rail Mounting Includes: 2 ea. Pivot Brackets (Left & Right), 1 ea. Tube (3/4 x 28"), 4 ea. Locknuts (1/4-20), and 4 ea. Screws, Hex Head (1/4-20 x 5/8").	1	1
NOTE: A - Includes: Items 2, Instruction Sheet, and Attaching Hardware. B - Not shown.					