


Bed Rails

Model IH5831 - 1/2 Length, Foot End

Model IH5831 - Discontinued 11/03



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	027040-NLA	NLA - Side Rail Assembly (1/2 Length) (Pair)	1	1
2		1113733-NLA	NLA - Side Rail, 1/2 Length - Right	1	1
2		1113734-NLA	NLA - Side Rail, 1/2 Length - Left	1	1
3		033786	Screw, Socket Head (#10-32 x 3/8")	1	4
4		027035	Tube, Side Rail	1	1
	B	020258	Wrench, Allen (5/32")	1	1
	B	026715-NLA	NLA - Instruction Sheet, Bed Rail	1	1
	B	026195	Hardware, Rail Mounting Includes: 2 ea. Pivot Bracket, 4 ea. Locknut (1/4-20), and 4 ea. Screw, Hex Head Cap (1/4-20 x 5/8").	1	1
NOTE: A - Includes: Items 2, Instruction sheet, and Attaching Hardware. B - Not Shown.					