


Accessories - Anti-Tip Assembly

Backrest Frame, Seat Frame
And Base Frame Shown
For Reference Only



Shown For
Reference Only

Accessories - Anti-Tip Assembly

Item Number	Note	Part Number	Description	Quantity Per			
				Package	Assembly		
1	A,C	1169817	Anti-Tip Assembly, Rear (Pair)	1	1		
2			NLA - Anti-Tipper	2	1		
4			NLA - Bracket, Anti-Tip - Right	1	1		
5			NLA - Bracket, Anti-Tip - Left	1	1		
7			NLA - Screw, Cap, Socket Head (1/4-28 x 1")	4	1		
8			NLA - Spacer, Bracket (1/4 x 1/2 x 7/16")	4	1		
1			B,C	1169818-NLA	NLA - Anti-Tip Assembly, Rear (Pair)	1	1
3					NLA - Anti-Tipper w/o Attaching Hardware	2	1
6	NLA - Bracket, Anti-Tipper, Black	2			1		
7	NLA - Screw, Cap, Socket Head (1/4-28 x 1")	4			1		
8			NLA - Spacer, Bracket (1/4 x 1/2 x 7/16")	4	1		
		See Chart	Anti-Tip Rear Seat To Floor Reference Chart	1	1		

NOTE: A - Includes items 2, 4, 5, 7, and 8.
 B - Includes items 3, 6, 7, and 8.
 C - See Chart titled "Anti-Tip Rear Seat to Floor Reference" for when to use which Anti-Tip Assembly.

Anti-Tip Rear Seat To Floor Reference Chart

Rear Seat to Floor	CRF22SP, CRF6722 (22" Wheel)	ZEROG24, CRF6724 SPLX24, FN1624 (24" Wheel)	ZEROG25, CRF6725 SPLX25, FN1625 ZERO26, CRF6726, FN1626 (25 & 26" Wheels)
14"	1169817	1169817	1169817
14-1/2"	1169817	1169817	1169817
15"	1169817	1169817	1169817
15-1/2"	1169817	1169817	1169817
16"	1169817	1169817	1169817
16-1/2"	1169818	1169817	1169817
17"	1169818	1169817	1169817
17-1/2"	1169818	1169818	1169817
18"	1169818	1169818	1169817
18-1/2"	1169818	1169818	1169818
19"	1169818	1169818	1169818
19-1/2"	1169818	1169818	1169818
20"	1169818	1169818	1169818