


Rollator (Four Wheeled)

Model 65800 - Sandstorm

Model 65800 (Sandstorm)
Denim Blue Finish



Rollator (Four Wheeled) Model 65800 - Sandstorm

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1076186	Grips, Hand	2	1
2	A	1076184-NLA	NLA - Kit, Hand Brake	1	2
3		1099372-NLA	NLA - Cover	1	1
4		1094528-NLA	NLA - Pad, Backrest	1	1
5	B	1099631-NLA	NLA - Kit, Backrest Tube, Denim Blue	1	1
6		1099367-NLA	NLA - Package, Anti-Rattle, Black (1")	2	1
7	C	1076175-NLA	NLA - Knob, Adjustable with Hardware	1	2
8		1099299-NLA	NLA - Seat	1	1
9	D	1099373-NLA	NLA - Kit, Seat Hardware	1	1
10	G	1110749-NLA	NLA - Kit, Connection Link w/ Hardware	1	2
11		1099369-NLA	NLA - Fork, Front	2	1
12	E	1099370-NLA	NLA - Kit, Wheel (5"), Front	2	1
13	F	1099371-NLA	NLA - Kit, Wheel - Rear (5")	2	1
14		1098988-NLA	NLA - Basket	1	1
15		1076174-NLA	NLA - Tray	1	1

NOTE: A - Includes Brake Handle, Brake Holder and Brake Cable Only.
B - Includes Back Tube and Snap Buttons.
C - Includes Knob and Bolt.
D - Includes Brackets and Screws.
E - Includes Caster, Bearing, Bolt, Washer and Nut.
F - Includes Caster, Bearing, Bolt, Washer, Spacer and Nut.
G - Includes Front and Rear Linkage, Hardware and Instruction Sheet.