


Drive Wheels for Free Wheel Hubs 4-Pole Motors Only



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	Item #'s 1, 2 and 3 Drive Wheels for Free Wheel Hubs	1	1
4		1104889	Rim, 8" w/ Hardware (14 x 3" Wheel)	1	1
4		1078246-NLA	NLA - Kit, 6" Rim w/ Hardware (14" x 4" Wheel)	1	1
	B,C	7930326-NLA	NLA - Package, Screw, Hex Head (1/4-20 x 3/4")	25	1
	B,C	7930333	Locknut (1/4-20)	25	1
	A,B	1024443-NLA	NLA - Screw, Hex Head (5/16-24 x 5/8")	1	4
	A,B	1024442-NLA	NLA - Washer, Lock (5/16 x 9/16 x 1/16")	1	4
	A,B	1024426-NLA	NLA - Nut, Hex (5/16-24)	1	4
NOTE: A - Included in rim with hardware part number 1078246 B - Not Shown C - Included in rim with hardware part number 1104889					

Item #'s 1, 2 and 3 Wheel Components

Wheel Size	Item #1 Wheel Complete Pneumatic	Item #1 Wheel Complete w/ Flat Free	Item #2 Tire	Item #2 Tire w/ Flat Free	Item #3 Tube
14" x 3"	1079298	1079301	1039652	1101754	1039651
14" x 4"	1079305-NLA	1079304	1042474	1042478	1042475