

Free Wheel Hub Kits (4 Pole Motors Before 8/17/10 US and Canada)

NOTE: Use kit numbers below to change Non-Free Wheel Hubs to Free Wheeling Hubs.

NOTE: All free wheel hub kits include items 2-4.

Free Wheel Hub Kit For (14" x 3" Wheel, Pneumatic)

Item Number	Part Number	Description	Quantity
1.	1079640-NLA See Note(s) A	Free Wheel Hub Kit 14" x 3" Wheel, Pneumatic	1
2.	1079004	Free Wheel Hub Assembly	2
3. See Note(s) B	1079295-NLA See Note(s) B	Instruction Sheet	1
4.	1081358	Wheel Assembly w/ Hub, 14" x 3" Pneumatic	2

NOTE: A - Kit includes items 2-4
B - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

Free Wheel Hub Kit For (14" x 3" Wheel w/ Flat Free)

Item Number	Part Number	Description	Quantity
1.	1079641-NLA See Note(s) A	Free Wheel Hub Kit 14" x 3" Wheel w/ Flat Free	1
2.	1079004	Free Wheel Hub Assembly	2
3. See Note(s) B	1079295-NLA	Instruction Sheet	1
4.	1081359	Wheel Assembly w/ Hub, 14" x 3" w/ Flat Free	2

NOTE: A - Kit includes items 2-4
B - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

Free Wheel Hub Kits (4 Pole Motors Before 8/17/10 US and Canada)

NOTE: Use kit numbers below to change Non-Free Wheel Hubs to Free Wheeling Hubs.

NOTE: All free wheel hub kits include items 2-4.

Free Wheel Hub Kit For (14" x 4" Wheel, Pneumatic)

Item Number	Part Number	Description	Quantity
1.	1079644-NLA See Note(s) A	Free Wheel Hub Kit 14" x 4" Wheel, Pneumatic	1
2.	1079004	Free Wheel Hub Assembly	2
3. See Note(s) B	1079295-NLA	Instruction Sheet	1
4.	1081362-NLA	Wheel Assembly w/ Hub, 14" x 4" Pneumatic	2

NOTE: A - Kit includes items 2-4
B - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

Free Wheel Hub Kit For (14" x 4" Wheel w/ Flat Free)

Item Number	Part Number	Description	Quantity
1.	1079645-NLA See Note(s) A	Free Wheel Hub Kit 14" x 4" Wheel, w/ Flat Free	1
2.	1079004	Free Wheel Hub Assembly	2
3. See Note(s) B	1079295-NLA	Instruction Sheet	1
4.	1081363	Wheel Assembly w/ Hub, 14" x 4" w/ Flat Free	2

NOTE: A - Kit includes items 2-4
B - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.