


Wheel Lock Assembly (4 Pole Motors After 8/16/10 US and Canada)



Wheel Lock Assembly (4 Pole Motors After 8/16/10 US and Canada)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1139885	Wheel Lock Assembly - Right	1	1
1	A	1139886	Wheel Lock Assembly - Left	1	1
2		1101725	Cap, Brake Lever End, Black	1	1
2		1130584	Package, Cap, Brake Lever End, Black	10	1
3	B	1169187	Kit, Wheel Lock Mounting Hardware	1	1
4			NLA - Screw, Hex Head w/ Patch (5/16-18 x 7/8")	1	2
5			Washer, Flat (5/16 x 5/8 x 1/16")	1	2
6			Locknut (5/16-18)	1	2

NOTE: A - Wheel Lock Assembly includes item #2 part number 1101725 Brake Lever End Cap.
B - Includes items 4-6 and accommodates one side only.