


Wheel Locks - Pull to Lock and Pull to Lock Extended



Wheel Locks - Pull to Lock and Pull to Lock Extended

| Item Number | Note | Part Number | Description | Quantity Per | |
|-----------------------------------|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1538079 | Wheel Lock Assembly w/ Hardware, Pull to Lock - Right | 1 | 1 |
| 1 | | 1538078 | Wheel Lock Assembly w/ Hardware, Pull to Lock - Left | 1 | 1 |
| 2 | | 1541957 | Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Right | 1 | 1 |
| 2 | | 1541956 | Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Left | 1 | 1 |
| 3 | | 1423107 | Kit, Hardware, Clamp Brake (Pair) | 1 | 1 |
| 4 | | | Screw (M6 x 1 x 35mm) | 1 | 4 |
| 5 | | | Mounting Clamp, Wheel Lock | 1 | 2 |
| NOTE: A - Includes items 4 and 5. | | | | | |