


Rollator (Four Wheeled)

Model 65900 - Stardust

Model 65900 (Stardust)
Jewel Green Finish



Rollator (Four Wheeled) Model 65900 - Stardust

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1076186	Grips, Hand	2	1
2	A	1076184-NLA	NLA - Kit, Hand Brake	1	2
3		1099372-NLA	NLA - Cover	1	1
4		1094528-NLA	NLA - Pad, Backrest	1	1
5		1099367-NLA	NLA - Package, Anti-Rattle, Black (1")	2	1
6	B	1076175-NLA	NLA - Knob, Adjustable with Hardware	1	2
7		1099299-NLA	NLA - Seat	1	1
8	C	1099373-NLA	NLA - Kit, Seat Hardware	1	1
9	E	1110749-NLA	NLA - Kit, Connection Link w/ Hardware	1	2
10		65406-NLA	NLA - Kit, Curb Climber	1	2
11	D	1076176-NLA	NLA - Kit, Wheel w/ Hardware	1	4
12		1082534-NLA	NLA - Cap, Plastic, Black	8	1
13		1094535-NLA	NLA - Fork, Front	1	2
14		1098988-NLA	NLA - Basket	1	1
15		1076174-NLA	NLA - Tray	1	1
16		1104754-NLA	NLA - Pad, Brake	1	2

NOTE: A - Includes Brake Handle, Brake Holder and Brake Cable Only.
B - Includes Knob and Bolt.
C - Includes Brackets, Plates and Screws.
D - Includes Caster, Bearings, Spacer, Nut, Screw and Caps.
E - Includes Front and Rear Linkage, Hardware and Instruction Sheet.