


Mid-Line Swing Bracket Assembly



Mid-Line Swing Bracket Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1090195-NLA	NLA - Mid-Line Joystick Assembly	1	1
2		1027274	Adjustment Lock Assembly, Black	1	3
3		1035528-NLA	NLA - Tube, Adjustment, Inside, Mid-Line, Black	1	1
4		1016778-NLA	NLA - Spring, Extension (9/16 x 1-3/8" L)	1	1
5		1027267-NLA	NLA - Half Clamp, Threaded Holes	1	1
6		1035860-NLA	NLA - Tube, Adjustment, Outside, Mid-Line, Black	1	1
7		1087070-NLA	NLA - Tube, Mid-Line Adjustable, Black	1	1
8		1027266	Half Clamp, Non-Threaded Holes	1	1
9		7930349-NLA	NLA - Screw, Hex Head (1/4-20 x 1-1/8")	25	1
10		1049582-NLA	NLA - Bracket, Joystick Mounting, Black	1	1
11		1027095	Plug Button, Black (7/8")	1	1

NOTE: A - Includes Items 2-11