


Casters - Front and Mounting Hardware



Casters - Front and Mounting Hardware

Item Number	Description	AT901 Composite Dark Grey	AT902 Composite Dark Grey	AT903 Composite	1292 Black Composite	1492 Composite Pneumatic
		Urethane (5 x 1")	Urethane (6 x 1")	Urethane (8 x 1-1/4")	Urethane (8 x 1-1/4")	w/ Flat Free (8 x1-1/4")
1	Wheel Assembly 1pc.	1029496 See Note(s) A	1034352 See Note(s) A, C	1058649 See Note(s) A	1009481 See Note(s) A	1030773
2	Tire	N/A	N/A	1038157	40176X040	40176X040
3	Tube/Insert	N/A	N/A	N/A	1010904-NLA	1010904-NLA
4	Bearing	N/A	N/A	1040102	1006049	1006049
5	Hub Spacer	N/A	N/A	1020592	1020592	1020592
6	Front Axle	N/A	7930362	57667P007	57667P007	57667P007
7	Axle Locknut	N/A	7930363-NLA	7930352-NLA	7930352-NLA	7930352-NLA
8	Axle Washer	N/A	N/A	N/A	7930312	7930312
	Axle Spacer See Note(s) B	N/A	1053202	N/A	N/A	N/A
NOTE: A - Includes tire, tube/insert, bearings, and spacer. B - Not Shown C - Consists of (2) Casters.						