

Casters - Front (Before 7/1/11)

Wheel Style	Wheel Size	Tire Type	Bearing	Bearing Spacer	Caster w/ Bearings	Tire	Tube/ Insert
AT900	3"	Urethane	40117X026	1026172	1032190-NLA	N/A	N/A
5ALUM	5"	Urethane	40117X026	1026172	1026952	N/A	N/A
AT901	5"	Urethane	40117X026	1026172	1029496	N/A	N/A
6ALUM	6"	Urethane	40117X026	1026172	1035206-NLA	N/A	N/A
AT902	6"	Urethane	N/A	N/A	1034352 See Note(s) A	N/A	N/A
6SKYWAY	6"	Urethane	40117X026	1026172	1026956	N/A	N/A
6PNEU	6"	Pneumatic	40117X026	1026172	1044040	1026010-NLA	1026011-NLA
AT903	8"	Urethane	1006049	1020592	1038161	1038157	N/A
AT1228	8 x 1-3/4"	Semi Pneumatic	1006049	1020592	1070811	N/A	N/A
AT1292	8 x 1-1/4"	Pneumatic	1006049	1020592	1009481	40176X040	40176X041
AT1292FF	8 x 1-1/4"	Pneumatic w/ Flat Free	1006049	1020592	1030773	40176X040	1010904-NLA

NOTE: A - Consists of (2) Casters.

6" Pneumatic Wheels with Flat-Free Inserts are No Longer Available (NLA). If a customer requires a pair, order a pair of 6" Composite Urethane Wheels. If a customer requires only one 6" Pneumatic Wheel with a Flat-Free insert, you will need to order a pair of 6" Composite Urethane Wheels so that both front wheels will match.