


Shroud Covers, Recline Actuator Curve Contoura w/ Cantilever Arms CGTR/CGTRE



Shroud Covers, Recline Actuator Curve Contoura w/ Cantilever Arms CGTR/CGTRE

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1169206-NLA	NLA - Shroud w/ Decal, Curve Contoura w/ Cantilever Arms	1	1
3	A	1169835-NLA	NLA - Kit, Curve Contoura Recline Actuator Shroud Mounting Hardware	1	1
4			Bracket, Shroud Rear Mount Clevis, Black	1	1
5			Screw, Socket Head w/ Patch (5/16-18 x 2")	1	2
6			Washer, Flat (5/16 x 5/8 x 1/16")	1	2
7			Locknut (5/16-18)	1	2
8			Screw, Thumb (#10-32 x 1/2")	1	3
NOTE: A - Includes items 4-8					
Item #2 Reflective Decal is not available for sale individually.					